



Forget everything you think you know about contact lenses...
Onefit is the game changer.

Onefit™ Scleral Lenses have been designed to provide crystal clear, and reliable vision performance with exceptional comfort. These lenses are an excellent alternative to soft or gas permeable (GP) lenses:

- Soft lenses can dry out causing discomfort, fluctuating vision and the lenses have the potential to pop out.
- GP lenses offer sharp vision, but may be less comfortable for some, may not be stable on the eye and might pop out easily.
- Scleral lenses offer the consistent, crisp vision of GP lenses with the hydrating comfort of soft lenses. They will not dry and are very unlikely to fall out of your eye!

Learn more at onefitlenses.com

Scleral lenses are a great option for just about anyone!


- Those who experience dryness, discomfort, or contact lens intolerance
- Athletes and active individuals – all eye shapes and conditions
- Those with Keratoconus, Astigmatism (toric), or Presbyopia (multifocal or bi-focal), and patients with less than optimum results from LASIK and Post-RK
- Normal corneas and even mild to highly irregular corneas

Benefits of wearing Onefit Scleral Lenses:

- Consistent all day comfort and hydration without risk of lens popping out during activity
- Crisp vision with no fluctuation, halos or glare
- Optimum oxygen to your cornea for long-term corneal health
- Lens and tear layer provide protection from wind, dust and debris

How do Onefit scleral lenses work?

Onefit lenses are made of a thin material that lets oxygen pass through the lens. They incorporate a unique design that eliminates irregularities of the shape on the front surface of your eye to correct your vision. Onefit scleral lenses are applied to the eye with a non-preserved saline solution that provides a thin cushion of fluid between the lens and your eye – providing a healthy and comfortable environment for your eye.



*Scleral lens vaults the cornea,
maintaining a cushion of tears*

"I have never been able to wear any type of contact lens because I have dry eyes and now I need bifocals to read. My doctor recently fit me with Onefit™ lenses and I love them! I wear them comfortably at least 15 hours a day, my vision is better than I see with my glasses, and I am able to read with them as well!"

S.B. – 52-year-old female, Atlanta, GA

Learn more at onefitlenses.com

